

Dializa

Nr	Përshkrimi i shkurtër	Gjendja aktuale ne stok
1	Low flux Biocompatible Membrane - High Performance, steam sterilization Effective Surface area (m2) 1.4 - deri 1.6; Ultrafiltration Coeff (ml/h - mmHg) me e madhe ose e barabart 10	
2	Low flux Biocompatible Membrane - High Performance, steam sterilization Effective Surface area (m2) 1.7 - deri 1.9; Ultrafiltration Coeff (ml/h - mmHg) me e madhe ose e barabart 12	
3	High flux Biocompatible Membrane - Haemodiafilter, steam sterilization Effective Surface area (m2) 1.4 - deri 1.6; Ultrafiltration Coeff (ml/h - mmHg) me e madhe ose e barabart 40	1080
4	High flux Biocompatible Membrane - Haemodiafilter, steam sterilization Effective Surface area (m2) 1.7 - deri 1.9; Ultrafiltration Coeff (ml/h - mmHg) me e madhe ose e barabart 50	5967
5	High flux Biocompatible Membrane - Haemodiafilter, steam sterilization Effective Surface area (m2) 2.0 - deri 2.2; Ultrafiltration Coeff (ml/h - mmHg) me e madhe ose e barabart 60	
6	ARTERIAL FISTULA NEEDLE 16 G 1.6mm Needle length 25 mm, tube length 150 mm, with rotating wing	45475
7	VENOUS FISTULA NEEDLE 16 G 1.6mm Needle length 25 mm, tube length 150 mm, with rotating wing	46010
8	TEMPORARY DOUBLE LUMEN CATHETER 11F, 200mm for vascular access in hemodialysis, suitable for all placement sites (vena subclavia, vena jugularis, vena femoralis)	200
9	TEMPORARY DOUBLE LUMEN CATHETER 11F, 150mm for vascular access in hemodialysis, suitable for all placement sites (vena subclavia, vena jugularis, vena femoralis)	5
10	DOUBLE LUMEN LONG TERM HEMODIALYSIS CATHETER FOR ADULT dime 14.5F x 23cm(nga Dializa)	67
11	KONCENTRAT I LENGHEM I BIKARBONATEVE 8.4% Kompatibil me aparaturën ekzistuese 4008S Kanister/Liter	0
12	KONCENTRAT I LENGHEM I BIKARBONATEVE 8.4% Kompatibil me aparaturën ekzistuese 4008S Kanister/Liter	
13	KONCENTRAT I LENGHEM ACIDIC ME PERMBAJTJE SI NE VIJIM: Na mmol/L = 138-140, K mmol/L = 2.00 , Ca mmol/L =1.25-1.50 , Acetate mmol/L = 3-6., Glucose-1gr/L Kanister/Liter	
14	KONCENTRAT I LENGHEM ACIDIK ME PERMBAJTJESI NE VIJIM: Na mmol/L = 138-140.00 , K mmol/L = 3.00 , Ca mmol/L =1.25-1.50 , Acetate mmol/L =3-6.00 , Glucose 0g/L Kanister/Liter	15620
15	NaCl (SALT) 25kg tab	20
16	SOLUCION PER DIALIZE PERITONEALE Solution for peritoneal dialysis,single use doble bag system made of non-PVC material for all system components, central control switch to regulate all treatment steps, Composition: 1.75 mmol Calcium and 1.5% glucose 2000ml	

17	SOLUCION PER DIALIZE PERITONEALE Solution for peritoneal dialysis, single use doble bag system made of non-PVC material for all system components, central control switch to regulate all treatment steps, Composition: 1.75 mmol Calcium and 4.25% glucose 2000ml	
18	SOLUCION PER DIALIZE PERITONEALE Solution for peritoneal dialysis, single use doble bag system made of non-PVC material for all system components, central control switch to regulate all treatment steps, Composition: 1.75 mmol Calcium and 2.3% glucose 2000ml	
19	- Filter i fluidit, kompatibel me aparatin (endotoxine retention capacity > 1.000.000, effective surface area 2.2 m2Kompatibil me apar.4008S	
20	- Filter i fluidit, kompatibel me aparatin (endotoxine retention capacity > 1.000.000, effective surface area 2.2 m2Kompatibil me apar. Dialog plus	
21	- Filter i fluidit, kompatibel me aparatin (endotoxine retention capacity > 1.000.000, effective surface area 2.2 m2Kompatibil me apar.AK 98	
22	Dezinfektues Ketrixh Cleancrt	
23	- Linjat A-V kompatible me aparturen. Diametri i dhomes venoze 18-22mm, Segmenti i pompës me diametër 8mm dhe gjatësi mbi 32 cm Kompatibil me aparaturat ekzistuese 4008S	9233
24	Tretësire dezinfektuese per aparat te hemodializes me baze te acidit citrik 20% paketim maksimal 5L (Dezinfektim I nxehtë)	
25	Tretësire dezinfektuese per aparat te hemodializes me baze te acidit citrik 50% paketim maksimal 5L (Dezinfektim I nxehtë)	
26	Tenckhof cathter double cuff kit	1
27	Transfer set	7
28	Catheter protection cap	810
29	CLAMP	30
30	PD solution heater plate	
31	Transfer set, multiuser adapter	5
32	Dezinfektues per duar	5
33	CAPD Scale	
34	SOLUCION PER DIALIZE PERITONEALE	2155
35	SOLUCION PER DIALIZE PERITONEALE	0
36	KONCENTRAT PLUHUR I BIKARBONATEVE 650GR BRAUN	1600
37	KONCENTRAT PLUHUR I BIKARBONATEVE 650GR FRESSENIUS	17965
38	SOLUCION PER DIALIZE PERITONEALE	341
39	DIACETRAL 4%	570
40	DIACETRAL 21%	0
41	DIACETRAL 50%	0
42	DIASEF FRESSENIUS	150
43	KONCENTRAT I LENGSEME ACIDIK ME KALIUM	0