

Dializa

Nr	Përshkrimi i shkurtër	Forma	Gjendja aktuale ne stok
1	Low flux Biocompatible Membrane - High Performance, steam sterilization Effective Surface area (m2) 1.4 - deri 1.6; Ultrafiltration Coeff (ml/h - mmHg) me e madhe ose e barabart 10	Cope	
2	Low flux Biocompatible Membrane - High Performance, steam sterilization Effective Surface area (m2) 1.7 - deri 1.9; Ultrafiltration Coeff (ml/h - mmHg) me e madhe ose e barabart 12	Cope	
3	High flux Biocompatible Membrane - Haemodiafilter, steam sterilization Effective Surface area (m2) 1.4 - deri 1.6; Ultrafiltration Coeff (ml/h - mmHg) me e madhe ose e barabart 40	Cope	1200
4	High flux Biocompatible Membrane - Haemodiafilter, steam sterilization Effective Surface area (m2) 1.7 - deri 1.9; Ultrafiltration Coeff (ml/h - mmHg) me e madhe ose e barabart 50	Cope	6887
5	High flux Biocompatible Membrane - Haemodiafilter, steam sterilization Effective Surface area (m2) 2.0 - deri 2.2; Ultrafiltration Coeff (ml/h - mmHg) me e madhe ose e barabart 60	Cope	
6	ARTERIAL FISTULA NEEDLE 16 G 1.6mm Needle length 25 mm, tube length 150 mm, with rotating wing	Cope	45805
7	VENOUS FISTULA NEEDLE 16 G 1.6mm Needle length 25 mm, tube length 150 mm, with rotating wing	Cope	46340
8	TEMPORARY DOUBLE LUMEN CATHETER 11F, 200mm for vascular access in hemodialysis, suitable for all placement sites (vena subclavia, vena jugularis, vena femoralis)	Cope	50
9	TEMPORARY DOUBLE LUMEN CATHETER 11F, 150mm for vascular access in hemodialysis, suitable for all placement sites (vena subclavia, vena jugularis, vena femoralis)	Cope	5
10	DOUBLE LUMEN LONG TERM HEMODIALYSIS CATHETER FOR ADULT dime 14.5F x 23cm(nga Dializa)	Cope	67
11	KONCENTRAT I LENGSHEM I BIKARBONATEVE 8.4% Kompatibil me aparaturën ekzistuese 4008S Kanister/Liter	Liter	0
12	KONCENTRAT I LENGSHEM I BIKARBONATEVE 8.4% Kompatibil me aparaturën ekzistuese 4008S Kanister/Liter	Liter	
13	KONCENTRAT I LENGSHEM ACIDIC ME PERMBAJTJE SI NE VIJIM: Na mmol/L = 138-140, K mmol/L = 2.00 , Ca mmol/L =1.25-1.50 , Acetate mmol/L = 3-6., Glucose-1gr/L Kanister/Liter	Liter	
14	KONCENTRAT I LENGSHEM ACIDIK ME PERMBAJTJESI NE VIJIM: Na mmol/L = 138-140.00 , K mmol/L = 3.00 , Ca mmol/L =1.25-1.50 , Acetate mmol/L =3-6.00 , Glucose 0g/L Kanister/Liter	Liter	23240
15	NaCl (SALT) 25kg tab	25kg Box	22
16	SOLUCION PER DIALIZE PERITONEALE Solution for peritoneal dialysis,single use doble bag system made of non-PVC material for all system components, central control switch to regulate all treatment steps, Composition: 1.75 mmol Calcium and 1.5% glucose 2000ml		
17	SOLUCION PER DIALIZE PERITONEALE Solution for peritoneal dialysis, single use doble bag system made of non-PVC material for all system components, central control switch to regulate all treatment steps, Composition: 1.75 mmol Calcium and 4.25% glucose 2000ml		
18	SOLUCION PER DIALIZE PERITONEALE Solution for peritoneal dialysis, single use doble bag system made of non-PVC material for all system components, central control switch to regulate all treatment steps, Composition: 1.75 mmol Calcium and 2.3% glucose 2000ml		
19	- Filter i fluidit, kompatibil me aparatin (endotoxine retention capacity > 1.000.000, effective surface area 2.2 m2Kompatibil me apar.4008S	Cope	
20	- Filter i fluidit, kompatibil me aparatin (endotoxine retention capacity > 1.000.000, effective surface area 2.2 m2Kompatibil me apar. Dialog plus	cope	

21	- Filter i fluidit, kompatibel me aparatin (endotoxine retention capacity > 1.000.000, effective surface area 2.2 m2Kompatibil me apar.AK 98	cope	
22	Dezinfektues Ketrihx Cleancrt		
23	- Linjat A-V compatibile me aparturen. Diametri i dhomezes venoze 18-22mm, Segmenti i pompës me diametër 8mm dhe gjatësi mbi 32 cm Kompatibil me aparaturat ekzistuese 4008S	Cope	9803
24	Tretësire dezinfektuese per aparat te hemodializes me baze te acidit citrik 20% paketim maksimal 5L (Dezinfektim I nxehtë)	Litra	
25	Tretësire dezinfektuese per aparat te hemodializes me baze te acidit citrik 50% paketim maksimal 5L (Dezinfektim I nxehtë)	Litra	
26	Tenckhof catheter double cuff kit		1
27	Transfer set		7
28	Catheter protection cap		810
29	CLAMP		30
30	PD solution heater plate		
31	Transfer set, multiuser adapter		5
32	Dezinfektues per duar		5
33	CAPD Scale		
34	SOLUCION PER DIALIZE PERITONEALE	1.36%	2155
35	SOLUCION PER DIALIZE PERITONEALE	2.27%	0
36	KONCENTRAT PLUHUR I BIKARBONATEVE 650GR BRAUN	QESE	400
37	KONCENTRAT PLUHUR I BIKARBONATEVE 650GR FRESENIUS	QESE	18085
38	SOLUCION PER DIALIZE PERITONEALE	3.86%	341
39	DIACETRAL 4%	L	570
40	DIACETRAL 21%	L	0
41	DIACETRAL 50%	L	60
42	DIASEF FRESENIUS	Cope	150
43	KONCENTRAT I LENGHEM ACIDIK ME KALIUM	LITER	0