

## GJENDJA E ANTIBIOTIKËVE

| NR | PERSHKRIMI                         | SASIA (Flakon x teste) |
|----|------------------------------------|------------------------|
| 1  | Amikacin 30 µg                     | 5x 50 teste            |
| 2  | Amphotericin B 20 µg               | 15x 50 teste           |
| 3  | Aztreonam 30 µg                    | 20 x 50 teste          |
| 4  | Azythromycin 15 µg                 | 28 x 50 teste          |
| 5  | Bacitracin 10 IU                   | 7 x 50 teste           |
| 6  | Cefalexin 30 µg                    | 5x 50 teste            |
| 7  | Cefepime 30 µg                     | 6 x 50 teste           |
| 8  | Cefixime 5 µg                      | 8 x 50 teste           |
| 9  | Cefotaxime 5 µg                    | 10 x 50 teste          |
| 10 | Cefoxitin 5 µg                     | 12 x 50 teste          |
| 11 | Ceftazidime 10µg                   | 8 x 50 teste           |
| 12 | Ceftriaxone 30µg                   | 12x 50 teste           |
| 13 | Cefuroxime 30 µg                   | 20x 50 teste           |
| 14 | Cloramphenicol 30µg                | 4 x 50 teste           |
| 15 | Ciprofloxacin 5µg                  | 10 x 50 teste          |
| 16 | Clindamycin 2µg                    | 17x 50 teste           |
| 17 | Colistin sulfate 10µg              | 4x 50 teste            |
| 18 | Ceftobiprole                       | 14x 50 teste           |
| 19 | Doripenem 10µg                     | 15 x 50 teste          |
| 20 | Eravacyclin 20µg                   | 6 x 50 teste           |
| 21 | Erythromycin 15µg                  | 26x 50 teste           |
| 22 | Fluconazole 25 µg                  | 17x 50 teste           |
| 23 | Flucytosine 10 µg                  | 17x 50 teste           |
| 24 | Fosfomicin 200µg                   | 17x 50 teste           |
| 25 | Fusidic acid 10µg                  | 21x 50 teste           |
| 26 | Gentamicin 30µg                    | 6x 50 teste            |
| 27 | Gentamicin 10 µg                   | 9x 50 teste            |
| 28 | Itraconazole 8µg                   | 18x 50 teste           |
| 29 | Ketonazole 10µg                    | 19 x 50 teste          |
| 30 | Levofloxacin 5µg                   | 12 x 50 teste          |
| 31 | Linezolid 10µg                     | 20x 50 teste           |
| 32 | Metronidazole 5µg                  | 6x 50 teste            |
| 33 | Miconazole 10µg                    | 16x 50 teste           |
| 34 | Moxifloxacin 5µg                   | 11x 50 teste           |
| 35 | Mupirocin 200µg                    | 12x 50 teste           |
| 36 | Nalidixic acid 30µg                | 10x 50 teste           |
| 37 | Nitrofurantoin 100µg               | 17x 50 teste           |
| 38 | Norfloxacin 10µg                   | 14x 50 teste           |
| 39 | Novobiocine 5µg                    | 10x 50 teste           |
| 40 | Nystatin 100 IU                    | 19 x 50 teste          |
| 41 | Optochin 5µg                       | 12 x 50 teste          |
| 42 | Oxacillin 1µg                      | 9 x 50 teste           |
| 43 | Penicillin G 1IU                   | 15 x 50 teste          |
| 44 | Piperacillin 30µg                  | 8 x 50 teste           |
| 45 | Piperacillin/Tazobactam 36 µg      | 9x 50 teste            |
| 46 | Quinopristin/Dalfopristin 15µg     | 14 x 50 teste          |
| 47 | Rifampicin 5µg                     | 21x 50 teste           |
| 48 | Teicoplanin 30µg                   | 9x 50 teste            |
| 49 | Tetracycline 30µg                  | 10 x 50 teste          |
| 50 | Tygecycline 15µg                   | 6x 50 teste            |
| 51 | Tobramycin 10µg                    | 6x 50 teste            |
| 52 | Trimethoprim/Sulfamethoxazole 25µg | 7x 50 teste            |
| 53 | Vancomycin 5µg                     | 22x 50 teste           |
| 54 | Voriconazole 1µg                   | 16x 50 teste           |
| 55 | Amoxicillin-clavulanic acid 30 µg  | 15x 50 teste           |
| 56 | Ampicillin 2 µg                    | 9 x 50 teste           |
| 57 | Ampicillin 10 µg                   | 13 x 50 teste          |
| 58 | Ampicillin-sulbactam 20 µg         | 11 x 50 teste          |
| 59 | Cefiderocol 30 µg                  | 15x 50 teste           |
| 60 | Caspofungin 5 µg                   | 10x 50 teste           |
| 61 | Ceftazidime/avibactam 14µg         | 5 x 50 teste           |

|                          |                                                                                                        |                       |
|--------------------------|--------------------------------------------------------------------------------------------------------|-----------------------|
| 62                       | Ceftolozane-Tazobactam 30 mcg                                                                          | 18 x 50 teste         |
| 63                       | Ertapenem 10µg                                                                                         | 12x 50 teste          |
| 64                       | Imipenem 10 µg                                                                                         | 5 x 50 teste          |
| 65                       | Meropenem 10µg                                                                                         | 4 x 50 teste          |
| 66                       | Meropenem/Vaborbactam 30 mcg                                                                           | 14x 50 teste          |
| 67                       | Colisitin Broth microdilution                                                                          | 45                    |
| 68                       | Compas colistin Broth                                                                                  | 16x4                  |
| 69                       | Compas Vancomycin Broth                                                                                | 2x8                   |
| 70                       | Compas antifungali Broth                                                                               | 18                    |
| 71                       | Coagulaze Test                                                                                         | 3                     |
| Nr                       | <b>Pershkrimi MIC TESTE</b>                                                                            | <b>SASIA ( PAKO )</b> |
| 1                        | Trimethoprim 0.016 - 256 mg/L                                                                          | 7x30                  |
| 2                        | Ampicillin 0.016-256 µg/mL                                                                             | 7x 30                 |
| 3                        | Cefiderocol 0.016-256 µg/mL                                                                            | 5 x 30                |
| 4                        | Ceftotaxime 0.002-32 µg/mL                                                                             | 1 x 30                |
| 5                        | Ceftazidime/avibactam 0.016-256 µg/mL                                                                  | 3 x 30                |
| 6                        | Ceftolozane/tazobactam 0.016-256 µg/mL                                                                 | 5x 30                 |
| 7                        | Ceftriaxone 0.002-32 µg/mL                                                                             | 2x 30                 |
| 8                        | Imipenem 0.002-32 µg/mL                                                                                | 9x 30                 |
| 9                        | Imipenem- relebactam 0.002/4 -32/4 mg/L                                                                | 4 x 30                |
| 10                       | Ketonazole 0.002-321 µg/mL                                                                             | 3x 30                 |
| 11                       | Meropenem/Vaborbactam 0.016-256 µg/mL                                                                  | 8x 30                 |
| 12                       | Penicillin 0.002-321 µg/mL                                                                             | 7x 30                 |
| 13                       | Vancomycin 0.016-256 µg/mL                                                                             | 4 x 30                |
| 14                       | Eravacycline 0.002-32 mg/L                                                                             | 4 x 30                |
| 15                       | Voriconazole 0.002-321 µg/mL                                                                           | 2 x 30                |
| 16                       | Tigecycline 0.016 - 256 mg/L                                                                           | 7x 30                 |
| 17                       | Meropenem 0.002-32 mg/L                                                                                | 8x 30                 |
| 18                       | Fluconazole 0.016-256 mg/L                                                                             | 3 x30                 |
| 19                       | Amphotericina B                                                                                        | 3X30                  |
| <b>STOKU VITEK</b>       |                                                                                                        |                       |
| Nr                       | <b>PERSHKRIMI</b>                                                                                      |                       |
| 1                        | 2GN for identification of enetobateriace                                                               | 195                   |
| 2                        | 2AST-N437 for suceptobility testing of clinicaly significant aerobic gr(-) bacili                      | 30                    |
| 3                        | 2AST-N437 for suceptobility testing of clinicaly significant aerobic gr(-) bacili                      | 120                   |
| 4                        | 2GP for identification of eneterococci, streptococci, staphylococci and a select gram (+) organisms    | 10                    |
| 5                        | AST-GP67 staphylococcus spp., enterococcus spp., streptococcus agalactia                               | 50                    |
| 6                        | 2 YST Identification and yeast like organisms                                                          | 42                    |
| 7                        | 2AST-YS 08 for susceptibilty testing of clinicaly significant yeast against speciedied antimicrobials. | 72                    |
| 8                        | ID Gr ANAEROBE                                                                                         | 19                    |
| 9                        | ID CAMP                                                                                                | 12                    |
|                          |                                                                                                        |                       |
|                          |                                                                                                        |                       |
|                          |                                                                                                        |                       |
| <b>STOKU PER VIDAS 3</b> |                                                                                                        |                       |
|                          | <b>TESTI</b>                                                                                           | <b>SASIA</b>          |
|                          | <b>A HCV</b>                                                                                           | <b>5</b>              |
| 1                        | Anti HEV IgG                                                                                           | 143                   |
| 2                        | HIV DUO                                                                                                | 410                   |
| 3                        | HBCT                                                                                                   | 325                   |
| 4                        | Anti HBS                                                                                               | 290                   |
| 5                        | HBCM                                                                                                   | 185                   |
| 6                        | HBE /HBET                                                                                              | 170                   |
| 7                        | TBRA                                                                                                   | 5                     |
| 8                        | HBsAg                                                                                                  | 350                   |
| 9                        | LYME IgG                                                                                               | 105                   |
| 10                       | LYME IgM                                                                                               | 104                   |
| 11                       | EBV VCAM IgM                                                                                           | 125                   |
| 12                       | EBNA IgG                                                                                               | 60                    |
| 13                       | CMV IgG                                                                                                | 243                   |
| 14                       | CMV IgM                                                                                                | 241                   |
| 15                       | RB IgG                                                                                                 | 262                   |
| 16                       | RB IgM                                                                                                 | 310                   |
| 17                       | TOX IgM                                                                                                | 301                   |
| 18                       | TOX IgG                                                                                                | 270                   |
| 19                       | HAVTotal                                                                                               | 67                    |

| STOKU PER HEMOKULTURA                       |                                        |          |
|---------------------------------------------|----------------------------------------|----------|
|                                             | TESTI                                  | SASIA    |
| 1                                           | HEMO AEROBE                            | 50       |
| 2                                           | HEMO AEROBE                            | 300      |
| 3                                           | HEMO PEDIATRIKE                        | 20       |
| 4                                           | HEMO PEDIATRIKE                        | 100      |
| 5                                           | HMOKULTUR ANAEROBE                     | 50       |
| 6                                           | HMOKULTUR ANAEROBE                     | 280      |
|                                             |                                        |          |
|                                             |                                        |          |
|                                             |                                        |          |
| STOKU PER TERRENE USHQVESE                  |                                        |          |
|                                             | PERSHKRIMI- TERRENET                   | SASIA    |
| 1                                           | Columbia Agar +5% blood                | 2200     |
| 2                                           | Chocolate Agar                         | 230      |
| 3                                           | Muller Hinton 2 Agar                   | 950      |
| 4                                           | Brucella Blood Agar                    | 6        |
| 5                                           | Chromogenic Candida                    | 300      |
| 6                                           | Mannitol Salt Agar                     | 25       |
| 7                                           | Mac Conkey Agar - Endo                 | 2600     |
| 8                                           | Yersinia Selective Agar                | 580      |
| 9                                           | SS Agar                                | 20       |
| 10                                          | XLD Agar                               | 450      |
| 11                                          | Mac Conkey Sorbitol Agar - SMACK       | 280      |
| 12                                          | Campylobacter CCDA Agar                | 120      |
| 13                                          | Thioglycollate Medium                  | 1300     |
| 14                                          | Chromogenic Sallmonella                | 10       |
| 15                                          | Selenite F Broth 10ml Tubes            | 1750     |
| 16                                          | Bile Aesculin Broth 10ml Tubes         | 1950     |
| 17                                          | Chromogenic Clostridium Difficile      | 100      |
| 18                                          | Chromogenic ESBL                       | 15       |
| 19                                          | Chromogenic MRSA                       | 10       |
| 20                                          | Chromogenic OXA-48                     | 30       |
| 21                                          | Chromogenic VRE                        | 18       |
| 22                                          | Sabouro Dextrose Agar                  | 25       |
| 23                                          | Mueller Hinton II Broth                | 1450     |
|                                             |                                        |          |
|                                             |                                        |          |
| STOKU PER TESTET RAPIDE IMUNOKROMATOGRAFIKE |                                        |          |
|                                             | TESTI                                  | SASIA    |
| 1                                           | ROTA+ADENOVIRUS                        | 1390     |
| 2                                           | CLOSTIRIUM TOX A+B                     | 60       |
| 3                                           | CRYPTO+GIARDIA+ENTOAMEBA               | 91       |
| 4                                           | NOROVIRUS G1 G2                        | 104      |
| 5                                           | ENTEROVIRUS                            | 65       |
| 6                                           | ASTROVIRUS                             | 80       |
| 7                                           | CAMPYLOBACTER                          | 98       |
| 8                                           | LEISHMANIA IGG/IGM                     | 28       |
| 9                                           | LEGIONELLA                             | 150      |
| 10                                          | MALARIA P.F/P.V/PAN                    | 35       |
| 11                                          | HELICO. PYLORI AG                      | 740      |
| 12                                          | GIARDIA                                | 95       |
| 13                                          | ENTAMOEBA                              | 28       |
| 14                                          | Leptospira                             | 10       |
| 15                                          | Cryptococcus                           | 39       |
| 16                                          | S.pneumoniae                           | 9        |
|                                             |                                        |          |
|                                             |                                        |          |
| STOKU I MATERIALIT SHPENZUES                |                                        |          |
|                                             |                                        | SASIA    |
| 1                                           | ANTI COLI I                            | 1X5 ML   |
| 2                                           | URINE SAMPLE CONTAINER STERILE 60      | 2x500    |
| 3                                           | URINE SAMPLE CONTAINER STERILE 100     | 2x500    |
| 4                                           | SAMPLE CONTAINER WITH SCREEN CAP 25 ML | 2X500    |
| 5                                           | ANAEROBIC W ZIP POUCH                  | 10X10    |
| 6                                           | LACTOPHENOL BLUE SOLUT.                | 2XML-1ML |
| 7                                           | ANTI COLI I                            | 1X5MML   |

|                                    |                                          |              |
|------------------------------------|------------------------------------------|--------------|
| 8                                  | ANTI COLI II                             | 1X1ML-1ML    |
| 9                                  | ANTICOLI III                             | 1X1ML-1ML    |
| 10                                 | ANTI COLI 0157                           | 1X1ML        |
| 11                                 | CRYO TUBES                               | 1X1ML-1ML    |
| 12                                 | ANTI SALMONELLA (a-e)                    | 3X1ML-1ML    |
| 13                                 | ANTI SALMONELLA (GR-D)                   | 4X1MML       |
| 14                                 | ANTI SALMONELLA(VI)                      | 2X1ML        |
| 15                                 | ANTI SALMONELLA (HGM)                    | 4X1MML       |
| 16                                 | ANTI SHIGELLA (POLY BOYDII 1-15)         | 2X1ML-1ML    |
| 17                                 | ANTI SHIGELLA (POLY DYSENTERIAE 1-10)    | 2X1ML-1ML    |
| 18                                 | ANTI SHIGELLA (POLY FLEXNERI DHE SONNEI) | 2X1ML-1ML    |
| 19                                 | ANTI SHIGELLA FLEXNERI                   | 2X1ML-1ML    |
| 20                                 | ANTI SHIGELLA SONNEI SANDE.              | 1X1ML-1ML    |
| 21                                 | THIOGLYCOLATE BROTH                      | 4X500        |
| 22                                 | SABOURAUD DEXTROSE AGAR                  | 3X500        |
| 23                                 | MAC CONKEY AGAR                          | 5X500        |
| 24                                 | MUELLER HINTON AGAR II                   | 6X500        |
| 25                                 | MUELLER HINTON AGAR                      | 5X500        |
| 26                                 | COAGULASE PLASMA REAGENT 5X4 ML          | 3X45         |
| 27                                 | OXIDASE TEST DISCS                       | 5            |
| 28                                 | GRAM STAIN CRYSTAL VIOLET                | 5X1L         |
| 29                                 | GRAM STAIN SAFRANINE                     | 4X1L         |
| 30                                 | GRAM STAIN IODINE LUGOL=SOLUTION         | 4x1L         |
| 31                                 | GRAM STAIN DECOLORIZER                   | 8X1L         |
| 32                                 | ANTI SALMONELLA I (A-E) + VI             | 1X1ML        |
| 33                                 | ANTI SALMONELLA Gr.D                     | 2X1ML        |
| 34                                 | ANTI SALMONELLA VI                       | 4XiML        |
| 35                                 | ANTI SALMONELLA Gr.H,G,M                 | 2X100        |
| 36                                 | COAGULASTE TEST                          | 12           |
| 37                                 | OPTOHINE ETST                            | 4X68         |
| 38                                 | COMPAS COLISTINE                         | 8            |
| 39                                 | AUTOCLAVE BAGS                           | 274X480      |
| 40                                 | PREVI COLOR RA2 A                        | 3            |
| 41                                 | PREVI COLOR RA2 B                        | 4            |
| 42                                 | PREVI COLOR RA2 C                        | 3            |
| 43                                 | KLIGER AGAR                              | 1X200        |
| 44                                 | UREA AGAR BASE                           | 1X300        |
| 45                                 | Inoculating lopss Polystrene 1m-10ml     | 12x1000      |
| 46                                 | Microscope slides frosted                | 150x50       |
| 47                                 | Microskop cover glasses                  | 25x100       |
| 48                                 | PETRI DISHES POLYSTYRENEC 90X17MM        | 5X500        |
| 49                                 | PETRI DISHES POLYSTYRENEC 60X15MM        | 4X500        |
| 50                                 | PETRI DISHES POLYSTERNEC 120X17          | 5X500        |
| 51                                 | CULTURE TUBES GLASS                      | 40X100       |
| 52                                 | SWAB STERILE (TRAMPON SWAP)              | 100X100      |
| 53                                 | TUBE 5MML POLYSTYRENE 12X75MM            | 10X500       |
| 54                                 | OXIDAS TEST                              | 10X30        |
| 55                                 | MICROAEROPHILIC BAGS 2-3LIT              | 14X10        |
| 56                                 | CRYO TUBES 2.0MML                        | 2X80         |
| 57                                 | AUTOCLAVES BOHAZART WASTE BAGGS 600X760  | 25X100       |
| 58                                 | AUTOCLAVES BOHAZART WASTE BAGGS 300X660  | 25X100       |
| 59                                 | STERILE PIPETTES 250 MIKROLITER          | 12X96        |
| 60                                 | STERIL PIPETTE TIPS 1000 MICROLITER      | 5X96         |
| 61                                 | ACRIDINE ORAMGE STAIN                    | 2X250        |
| 62                                 | IMMERSION OIL                            | 6X100        |
| 63                                 | Inoculating lopss Polystrene 10ml        | 30x1000      |
| <b>GJENDJA E STOKUT PER AMPLEX</b> |                                          |              |
|                                    |                                          | <b>SASIA</b> |
|                                    | EAZYLPEX SUPER BUG ACINETO               | 10           |
|                                    | EAZYPLEX SUPER BUG MCR-I                 | 12           |
|                                    | EAZYLPEX SUPER BUG COMPLET A             | 3            |
|                                    | EAZYPLEX CSF DIRECT V                    | 1            |
|                                    | EAZYPLEX CSF DIRECT M                    | 1            |
|                                    | KLEBSIELLA PNEUMONIAE                    | 21           |
|                                    | EAZYPLEX SUPER COMPLET C                 | 19           |
|                                    |                                          |              |

| GJENDJA E STOKUT PER BCU CURETIS |                                  |       |
|----------------------------------|----------------------------------|-------|
|                                  |                                  | SASIA |
|                                  | BCU BLOOD CULTURE CATRIDGE SET   | 15    |
|                                  | HPN PNEUMINIAE                   | 4     |
|                                  | MASTER MIX TUBES                 | 28    |
|                                  | <b>GEN EXPERT</b>                |       |
|                                  | HBV-DNA                          | 49    |
|                                  | HCV-RNA                          | 18    |
|                                  |                                  |       |
|                                  |                                  |       |
|                                  |                                  |       |
|                                  | PNEUMOSLIDE                      |       |
|                                  | PNEUMISLIDE IGM                  | 29    |
|                                  | PNEUMISLIDE IGG                  | 37    |
|                                  |                                  |       |
|                                  |                                  |       |
|                                  |                                  |       |
|                                  | PHOENIX                          | Sasia |
|                                  | Panel phoenix NMIC 502           | 5     |
|                                  | Panel phoenix PMIC 96            | 49    |
|                                  | Panel phoenix NMIC 402           | 42    |
|                                  | Panel phoenix PMIC 90            | 50    |
|                                  | Tube Phoenix AST broth           | 162   |
|                                  | Phoenix AST indicator            | 26    |
|                                  | AST broth                        | 185   |
|                                  | Phoenix Pipette tips with filter | 3X163 |
|                                  | Tube phoenix inculum Broth       | 151   |
|                                  |                                  |       |
|                                  |                                  |       |

|                        |       |
|------------------------|-------|
| DIASORIN               | Sasia |
| TOXO IGM               | 200   |
| TOXO IgG II            | 100   |
| RUBELLA IgG II         | 100   |
| RUBELLA IgM            | 200   |
| CMV IgG II             | 200   |
| CMV IgM II             | 200   |
| HSV -1/2 IgG           | 200   |
| HSV-1/2 IgM            | 200   |
| HSV-2 IgG              | 200   |
| HSV-1 IgG              | 200   |
| Measles IgG            | 100   |
| Measles IgM            | 50    |
| Control Toxo IgM       | 4     |
| Control Toxo IgG II    | 4     |
| Control Rubella IgG II | 4     |
| Control Rubella IgM    | 4     |
| Control CMV IgG II     | 4     |
| Control CMV IgM II     | 4     |
| Control HSV-1/2 IgG    | 4     |
| Control HSV-1/2 IgM    | 4     |
| Control HSV-2 IgG      | 4     |
| Control HSV-1 IgG      | 4     |
| Control Measles IgG    | 4     |
| Control Measles IgM    | 4     |
| Cleaning Kit           | 12    |
| Wash system liquid     | 120   |
| Starter Kit            | 6000  |
| DiTi                   | 13824 |
| Cuvettes               | 14400 |
| Solid Waste Beg        | 100   |